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TM100 Restraint Lateral Splints Trapezometacarpal Covered The Detention Lateral Splints TM150 trapezometacarpal Indicated in cases that require a functional immobilization of the thumb. Suitable stage joint pain, whether caused by inflammation, trauma, after surgery or noninvasive treatments below. Use the deviations caused by delay rheumatic deformities.

PLACEMENT: Remove the tape around the ferrule, insert the thumb hole distally, attaching the brace so that it covers the thenar. Set the tape in position and length to the wrist. The joint motion should be trapezometacarpal is limited. Flexion of the interphalangeal joint of the thumb should not be restricted. If so, see your prosthetist will proceed to cut the brace properly.

SIZE (WRIST CONTOUR CM.) S (<17), M (17-19) L (> 19)
SIDES: right and left.

De Quervain Ferula DQ100 De Quervain Ferula DQ150 Covered Indicated in cases that require immobilization of the joint preventing trapezometacarpal flexo-extensor movement of the wrist joint. Maintains the thumb abducted, in its functional position. Indicated in painful and inflammatory processes that require immobilization of these joints, stenosing tenosynovitis of the abductors of the thumb, rheumatism, associated with anti-inflammatory treatments, such as postoperative interventions performed at this level.

PLACEMENT: Remove the tape around the ferrule, insert the thumb hole distally, attaching the brace so that it covers the thenar. Adjust the straps in position and length in the wrist and forearm. Make sure the wrist ulnar movement, and movement of the joint trapezometacarpal are limited. Flexion of the interphalangeal joint of the thumb should not be restricted. If so, see your prosthetist will proceed to cut the brace properly.

SIZE (WRIST CONTOUR / LENGTH IN CM.) S (<17/19), M (17-19/21) L (> 19/23)
SIDES: right and left.

MM100 Flat Hand Splint Postural Indicated in cases that require a neutral position of the wrist and fingers. Cases of hemiplegia, spastic set of wrist and finger joints, bedridden for long periods in cases of acute rheumatoid arthritis pain relief and align joints in post-operative period in burned to minimize scar formation, and rehabilitation of the extension of the MCP

and IF, etc.
PLACEMENT: Remove the tape around the splint placed on the ventral forearm in the body of the splint with the wrist in neutral position and fingers resting on the flat, the film develops between the palmar finger and thumb, setting the fingers, is inserted into the slot on the side and crosses the area near the wrist, securing the palm. Adjust the straps in position and length in the palm, wrist and forearm.

SIZE (WIDTH IN CM PALMAR.) S (6-7), M (7-8.5) L (8.5-10)
SIDES: right and left.

MM200 Functional Hand Splint Postural Indicated in cases that require a neutral position of the wrist and fingers functional thumb opposition. Cases of hemiplegia, bedridden for long periods in cases of acute rheumatoid arthritis pain relief and align joints in post-operative period to minimize the loss of functional mobility. **PLACEMENT:** Remove the tape around the splint placed on the ventral forearm in the body of the splint with the wrist in neutral position and fingers in a functional position on the plastic zone and the thumb in opposition, as suggested splint, adjusting the tape at that level. Adjust the straps in position and length in the wrist and forearm.

SIZE (width in cm palmar.) S (6-7), M (7-8.5) L (8.5-10)
SIDES: right and left.

MM500 Wrist Splint Immobiliser Indicated in cases that require a secure immobilization of the wrist joint. Trauma cases for wrist-joint without fracture, after the period of detention for non-surgical plaster in fractures, injuries to the distal radio, post-operative synovitis and synovial cysts of the wrist, if you remove the external fixator the Colles-Pouteau fracture, osteoarthritis, degenerative and recurrent sprains of the wrist... **PLACEMENT:** Remove the tape around the splint placed on the ventral forearm in the body of the splint with the wrist in neutral position, inserting the thumb through the side. Adjust the straps in position and length in the palm, wrist and forearm, ensuring that the flap remains inside the body of the splint.

SIZE (wrist contour cm.) S (14-17), M (17-20) L (20-23)
SIDES: right and left.

MM600 Thomas Splint With Thumb Abduction Indicated as active splint in the treatment of the hand affected by radial nerve palsy

and rehabilitation of hand trauma or post surgery.
PLACEMENT: Adapting the dorsal splint for complying with the tapes, directly or plaster. The length of the stretch determines the strength of action, both in extension of the fingers and in thumb abduction.

ONE SIZE: Conformable and adjustable straps
SIDES: right and left.

ADAPTATION: The wider tape preferably be placed in the wrist to improve fit, but can be exchanged if necessary. The tapes allow it to cut to fit the patient to avoid excessive volume of the same, and can be removed for cleaning them. In the event that the brace is not perfectly coupled to the morphology of the hand, the prosthetist can proceed with cutting edges and / or its conformation by applying hot air.

Precautions And Hygiene It should prevent the stiffening of the joints from prolonged use, by prescription must be combined use with mobility exercises, postural changes, etc. Do not exert excessive pressure on the settings of the tapes, mainly in wrist and forearm. Do not expose the splint to a flame or heat source at high temperatures, it may deform or damage it. Washable with mild soap and water. Air dry, or towel.

Precautions Follow the manufacturer's directions for proper orthosis placement. Follow the instructions or recommendations of the prescriber or healthcare personnel. For a personalized adaptation, consult a specialist. If there is ongoing pain, it is advisable to contact a healthcare professional. If an allergic reaction should occur remove the product and consult with the prescribing physician. In case of minor annoyances caused by excessive sweating we recommend the use of a tubular cotton. Positioned such that compression is not excessive. Do not expose the garment to a flame or heat source at very high temperatures. Store at room temperature, preferably in the original packaging. Guarantee within the framework of laws.

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en

Regulations The medical device complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices and the respective national regulations. Class I medical device, non-sterile, without measurement function.

Notice To dispose of the product, use the solid waste container. For the disposal of the container, depending on its material, use the paper / cardboard container or the plastic container. In any case, you must strictly comply with the regulations of the country. Any serious incident related to the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and / or the patient are established.

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